

Human Adiponectin

Mouse Monoclonal Antibody, clone: Ha6

Catalog Number: 21010

Size: 100 µg
Host: Mouse
Isotype: Mouse IgG
Immunogen: HEK-293 derived recombinant human adiponectin (Cat. No. 41013)
Species reactivity: Human. Not yet tested in other species
Purification method: Affinity chromatography on a column with immobilized protein G.
Formulation: Solution in PBS
Storage: -20°C

Application:

This antibody can be used as a capture antibody in a human adiponectin ELISA in combination with monoclonal anti-human adiponectin antibody (Cat. No.: 21011)

Related products:

Cat. No.	Description	Host	Reactivity	Applications
21010	Adiponectin, MAb, clone Ha6	Mouse	Human	ELISA capture
21011	Adiponectin, MAb, clone Ha7	Mouse	Human	ELISA detection
41013	Recombinant human adiponectin, FLAG-tagged (Mammalian)			
41014	Recombinant human adiponectin, His-tagged (E.coli)			
31010	Human adiponectin ELISA kit		1.55-100 ng/ml	
31011	Human adiponectin ELISA kit, one-step		1.55-100 ng/ml	