

## Human Adiponectin

### Mouse Monoclonal Antibody, clone: Ha7

**Catalog Number:** 21011

**Size:** 100 µg  
**Host:** Mouse  
**Isotype:** Mouse IgG  
**Immunogen:** HEK-293 derived recombinant human adiponectin (Cat. No. 41013)  
**Species reactivity:** Human. Not yet tested in other species  
**Purification method:** Affinity chromatography on a column with immobilized protein G.  
**Formulation:** Solution in PBS  
**Storage:** -20°C

**Application:**

This antibody can be used as a detection antibody in a human adiponectin ELISA in combination with monoclonal anti-human adiponectin antibody (Cat. No.: 21010)

**Related products:**

Cat. No.	Description	Host	Reactivity	Applications
21011	Adiponectin, MAb, clone Ha7	Mouse	Human	ELISA detection
21010	Adiponectin, MAb, clone Ha6	Mouse	Human	ELISA capture
41013	Recombinant human adiponectin, FLAG-tagged (Mammalian)			
41014	Recombinant human adiponectin, His-tagged (E.coli)			
31010	Human adiponectin ELISA kit		1.55-100 ng/ml	
31011	Human adiponectin ELISA kit, one-step		1.55-100 ng/ml	