

## Human Lipocalin-2

### Mouse Monoclonal Antibody, clone: 2G1

**Catalog Number:** 21050

**Size:** 100 µg

**Isotype:** Mouse IgG1

**Immunogen:** *E.coli*-derived recombinant human lipocalin-2 (Cat. No: 41050).

**Species reactivity:** Human.

**Purification method:** Affinity chromatography on a column with immobilized protein G.

**Formulation:** Solution in PBS.

**Storage:** -20°C

**Application:**

This antibody can be used as a capture antibody in a human Lipocalin-2 ELISA in combination with monoclonal anti-human Lipocalin-2 antibody clone 2C3 (Cat. No: 21051).

**Related products:**

Cat. No	Description
41050	Recombinant human lipocalin-2
21051	Human lipocalin-2, MAb, clone 2C3
31050	Human lipocalin-2 Elisa Kit