

Human Lipocalin-2

HRP conjugated mouse monoclonal antibody, clone: 2C3

Catalog Number: 21051-HRP

Size: 100 µg
Host: Mouse
Isotype: Mouse IgG1
Immunogen: E.coli-derived recombinant human lipocalin-2 (Cat. No. 41050)
Species reactivity: Human. Not yet tested in other species
Purification method: Affinity chromatography on a column with immobilized protein G.
Formulation: Solution in PBS.
Storage: -20°C

Application:

This antibody can be used as a detection antibody in a human Lipocalin-2 ELISA in combination with monoclonal anti-human Lipocalin-2 antibody (Cat. No.: 21050)

Related products:

Cat. No	Description
41050	Recombinant human lipocalin-2
21050	Human lipocalin-2, MAb, clone 2G1
31050	Human lipocalin-2 Elisa Kit