

Human Retinol Binding Protein 4 (RBP4)

Mouse Monoclonal Antibody, clone: 5H9

Catalog Number: 21060

Size: 100 µg
Host: Mouse
Isotype: Mouse IgG
Immunogen: E. coli derived recombinant human RBP4 (Cat. No. 41060)
Species reactivity: Human. Not yet tested in other species
Purification method: Affinity chromatography on a column with immobilized protein G.
Formulation: Solution in PBS
Storage: -20°C

Application:

This antibody can be used as a capture antibody in a human retinol binding protein 4 ELISA in combination with monoclonal anti-human retinol binding protein 4 antibody (Cat. No.: 21061)

Related products:

Cat. No.	Description	Host	Reactivity	Applications
21060	Retinol binding protein 4, MAb, clone 5H9	Mouse	Human	ELISA capture
21061	Retinol binding protein 4, MAb, clone 1A2	Mouse	Human	ELISA detection
41060	Recombinant human retinol binding protein 4 (E.coli)			
31060	Human retinol binding protein 4 ELISA kit		6.25-200 ng/ml	