

## Human Plasminogen Activator Inhibitor type-1 (PAI-1)

### Mouse Monoclonal Antibody, clone: 3H3

**Catalog Number:** 21070

**Size:** 100 µg  
**Host:** Mouse  
**Isotype:** Mouse IgG  
**Immunogen:** E. coli derived recombinant human PAI-1 (Cat. No. 41070)  
**Species reactivity:** Human. Not yet tested in other species  
**Purification method:** Affinity chromatography on a column with immobilized protein G.  
**Formulation:** Solution in PBS  
**Storage:** -20°C

**Application:**

This antibody can be used as a capture antibody in a human PAI-1 ELISA in combination with polyclonal anti-human PAI-1 antibody (Cat. No.: 11070)

**Related products:**

Cat. No.	Description	Host	Reactivity	Applications
21070	PAI-1, MAb, clone 3H3	Mouse	Human	ELISA capture
11070	PAI-1, PAb	Rabbit	Human	ELISA, WB
41070	Recombinant human PAI-1 (E.coli)			
31070	Human PAI-1 ELISA kit		0.031-2 ng/ml	