

Polyclonal Antibody against Mouse Retinol protein 4(RBP4)

Catalog Number: 12060

Size: 100 ug

Host: Rabbit

Introduction to the Molecule

Retinol binding protein 4(RBP4), originally known as a specific transport of retinol in blood, is also a novel inflammatory and insulin resistance marker. Serum RBP4 levels are elevated in insulin resistant mice and humans with obesity and type 2 diabetes. Animal experiments found that increased secretion of RBP4 might reduce insulin-dependent glucose uptake by muscle tissue by reducing the activity of PI(3)K (phosphoinositide 3-kinase), and increased hepatic glucose output by increasing the expression of the enzyme PEPCK2. Studies suggested that elevated serum RBP4 was associated with components of the metabolic syndrome, including increased body-mass index, waist-to-hip ratio, serum triglyceride levels, and systolic blood pressure and decreased high-density lipoprotein cholesterol levels. Furthermore, circulating RBP4 concentrations were associated with subclinical cardiovascular disease, which imply that RBP4 could be involved in the development of atherosclerosis.

Isotype/Preparation:

Rabbit specific IgG was purified by Mouse RBP4 affinity chromatography

Immunogen:

Recombinant full-length mouse RBP4 expressed in E.coli.

Specificity:

The antibody detects mouse RBP4.

Formulation:

Solution in PBS. Store at -20°C . For long-term storage, aliquot and freeze at -70°C . Avoid repeated freeze/defrost cycles.

Application/Usage:

Western blot - This antibody can be used as primary antibody in western blot assay to detect mouse RBP4

ELISA - This antibody can be used in ELISA with the appropriate secondary reagents to detect mouse RBP4

Quality Control Test

BCA to determine quantity of the antibody

Contact Us

- Website: www.torontobioscience.com
- E-mail: sales@torontobioscience.com